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
【Storage of hydrogen gas(brief version)】

If the electrolyte solution is acidic ($\text{pH} < 2$), it can be used as inner solution of RHE directly. Otherwise, the 013375 Double junction holder should be used together with the RHE to maintain the proper potential of RHE.

You can check more detail manual in our website.

URL: usr.bas.jp/dl_sub/?id=529d7ecacf927

1. Connect the Tubing clamp(b), Female lure fitting(c), Disposable syringe(d) and Silicone tube(e) to RHE reversible hydrogen electrode(a) like left figure.

2. Insert the RHE(a) and a Pt counter electrode(options; 012961, 002222 etc.) in acidic solution.

3. Leave the tubing clamp(b) at open state. Suck the acidic solution gently by using the disposable syringe(d) until the solution reaches near the female lure fitting(c) then close the tubing clamp(b) firmly.

4. Connect the working electrode cable of your potentiostat to RHE, and connect the reference and counter electrode cable to the Pt counter electrode. Electrolyze the acidic solution until 2 coulomb electric charges pass. (For example: apply -3V vs. ref to WE for about 3 minutes in 1.2 M HCl solution.)

【CAUTION】

- This reference electrode should be used in aqueous solution. Do not use it in organic solvent
- The electrolysis current might be over 10mA. Please check the Max. current range of the potentiostat before electrolysis.
- To avoid breaking the glass body, the electrode should be protected from strong shock.
- When you will use with RRDE-3A, please check the position of the RHE. The RHE might be broken by collision with the motor unit of the RRDE-3A.
- Use the electrode at room temperature and atmospheric pressure.
- After use, rinse well inside of RHE by pure water.
- We do not guarantee the product, if you applied the electrode against the above notes.